


#### **Cheat Sheet**

# **Key Features**

- Single Point Adjustment
- Dental Grade Cr-Co Alloy (Ni Free)
- Snap-Fit<sup>™</sup> or Perfect-Fit<sup>™</sup> Trays
- 15mm protrusion range
- 10mm lateral excursion range
- Maximum lingual tongue space
- Easy to Prescribe, Deliver & Adjust


### **Protrusion Range**

• 15 mm w/ three hooks



## **Adjusting Protrusion**

- Advance counter-clockwise, facing patient
- Half turn = 0.25 mm
- More than 5mm, change hook




## **Key Benefits**

- Proven effective for 20 years with over 35 peer reviewed studies
- Preferred & trusted by the US Government
- Medicare approved and FDA cleared
- Add TAP PAP CS for Combo Therapy
- Low Chair Time:
  - o Fits 1st time & Patient can adjust

#### **Hardware Options**

- Quick Release Bar
- Medicare Bar
- 4-10 mm Lateral Excursion



## **Tray Options**

- TL Snap-Fit
- ThermAcryl Perfect-Fit
- 5-15mm Vertical height
- Posterior Stops
- Heavy Bruxer Option
- Compliance Chip add-on

## Standard Design

- Quick release bar
- Set at 10mm forward & 5mm backward
- About 10mm lateral excursion
- About 5mm of vertical in the anterior
- TL Snap-Fit Tray

comparison chart	dream TAP	Herbst	Dorsal Fin	Tray Style	
Single Point Adjustment	Yes	No	No	No	Bi-lateral adjustment can cause TMD problems
Maximum Protrusion Range	15 mm	5-7 mm	5-8.5 mm	5-15 mm	Bigger range makes it easy to hit the protrusion target
Micro-Adjustments < .5 mm increments	Yes	Yes	Yes	No	Fine tuning makes it easy to hit the protrusion target
Lateral Excursion Range	4-10 mm	2-4 mm	0 mm	3-4 mm	Greater lateral excursion is more comfortable for patient
Hardware Material	Cr-Co Ni & Be Free	SS	SS	Polymer	Dental grade Cr-Co is the best alloy for the oral environment, bar none
Tray Material	Snap-Fit Polycarbonate or Perfect-Fit ThermAcryl	Acrylic w/ball clasps, flex liner optional	Acrylic w/ball clasps, flex liner optional	Thermoplastic or Laser Sintered Polyamide	Flexible trays create greater retention without risking teeth-moving ball clasps. ThermAcryl allows for a perfect fit in the chair
Prevents Mouth Breathing	Yes	With hooks & rubber bands	With hooks & rubber bands	No	Eliminating mouth breathing improves efficacy of any appliance
Compliance Chip Available	Yes	Yes	Yes	Yes	Future DOT requirement
Professional pre-fab device available	Yes	No	No	No	myTAP® can be dispensed in the practice in 15min.
Free Bite repositioner	Yes	No	No	No	AM Aligner® included in every case
Add-on CPAP Mask for Combo Therapy	Yes	No	No	No	Add a mask to treat severe cases. Also improves MD relationships
Peer Reviewed Studies	36	6	4	4	TAP Appliances are the most studied, by far
ISO 13485 Certified Manufacturing	Yes	No	No*	No	* Dynaflex is ISO 13485 Certified
Supplier to the US Armed Services & VA	Yes	No	No	No	Exclusive appliance for the Army, Navy and Air Force
Number of Supplying Labs	150	12	12	150	Use your favorite lab to build your dreamTAPs
PDAC Approved for Medicare	Yes	Yes	No	No	Medicare is about 30-40% of the market. This is a must

**Dorsal Fin:** Somnodent, Respire, ClearDream, Moses and many more....

Tray Style: EMA, Norval, Pantera, Silent Night